

## **Wastewater - Al-Fe-Mn-Mo-Pb-Sn-Zn: ppm-level**

Art. ID LM-VKI-QC-METAL-HL1  
Unit 1 ampoule  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2031 II E2

### Description

0,25L prepared reference material per ampoule

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	2	mg/L		
	Iron (Fe)	[7439-89-6]	3	mg/L		
	Manganese (Mn)	[7439-96-5]	2	mg/L		
	Molybdenum (Mo)	[7439-98-7]	10	mg/L		
	Lead (Pb)	[7439-92-1]	10	mg/L		
	Tin (Sn)	[7440-31-5]	10	mg/L		
	Zinc (Zn)	[7440-66-6]	0,5	mg/L		