

Wastewater - Ag-Ba-Cd-Co-Cr-Cu-Ni-Sr: ppm-level

Art. ID LM-VKI-QC-METAL-HL2
Unit 1 ampoule
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2031 II E2

Description

0,25L prepared reference material per ampoule

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	2	mg/L		
	Barium (Ba)	[7440-39-3]	2	mg/L		
	Cadmium (Cd)	[7440-43-9]	1	mg/L		
	Cobalt (Co)	[7440-48-4]	0,5	mg/L		
	Chromium (Cr)	[7440-47-3]	4	mg/L		
	Copper (Cu)	[7440-50-8]	4	mg/L		
	Nickel (Ni)	[7440-02-0]	2	mg/L		
	Strontium (Sr)	[7440-24-6]	5	mg/L		