

Fluorinated compounds in Methanol

Art. ID	LM24-SL-70025
Unit	5 x 5 mL
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation 24 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Perfluoroheptanoic acid (PFHpA)	[375-85-9]	20	µg/mL		
	Perfluoro-n-octanoic ac id	[335-67-1]	20	µg/mL		
	Perfluoro-n-nonanoic ac id	[375-95-1]	20	µg/mL		
Perfluorobutanesulfonic Acid	Nonafluorobutane-1-sulf onic acid	[375-73-5]	20	µg/mL		
	Perfluoro-n-octane sulf onic acid	[1763-23-1]	20	µg/mL		
	Perfluorooctane sulfona mide (PFOSA)	[754-91-6]	20	µg/mL		
	Perfluorohexanoic acid (PFHxA)	[307-24-4]	20	µg/mL		
	Perfluorodecanoic acid (PFDA)	[335-76-2]	20	µg/mL		
2-(Perfluorohexyl)ethan e-1-sulfonic Acid	1H,1H,2H,2H-Perfluoro-n -octane sulfonic acid	[27619-97-2]	20	µg/mL		
	Perfluorobutanoic acid (PFBA)	[375-22-4]	20	µg/mL		
	Perfluoro-n-pentanoic a cid	[2706-90-3]	20	µg/mL		
	Tridecafluorohexane-1-s ulfonic acid potassium salt (Perfluoro-n-hexan e sulfonic acid potassi um salt)		20	µg/mL		