

## Water QC Reference Material - NO3/NH4/PO4

Art. ID	IV-QCP-NUT-1-20ML
Unit	20 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

Lot Number: K2-NUT01126

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified Value based on diluted solution (1:200)	Ammonium (NH <sub>4</sub> <sup>+</sup> ) as N		9,91 ± 0,02	mg/L	Ion Chromatography	
Certified Value based on diluted solution (1:200)	Nitrate (NO <sub>3</sub> <sup>-</sup> ) as N		2,992 ± 0,011	mg/L	Ion Chromatography	
Certified Value based on diluted solution (1:200)	Orthophosphate as P		5,003 ± 0,011	mg/L	Ion Chromatography	
Certified density of solution at 20 ± 1 °C	Density		1,006	g/mL		