

Water QC Reference Material - NO₃/NH₄/PO₄

Art. ID IV-QCP-NUT-1-20ML
Unit 20 mL
Deliverydetails No Dangerous Good /not restricted

Description

Lot Number: K2-NUT01126

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified Value based on diluted solution (1:200)	Ammonium (NH ₄ ⁺) as N		9,91 ± 0,02	mg/L	Ion Chromatography	
Certified Value based on diluted solution (1:200)	Nitrate (NO ₃ ⁻) as N		2,992 ± 0,011	mg/L	Ion Chromatography	
Certified Value based on diluted solution (1:200)	Orthophosphate as P		5,003 ± 0,011	mg/L	Ion Chromatography	
Certified density of solution at 20 ± 1 °C	Density		1,006	g/mL		