

Electrolytic Copper - chips

Art. ID ERM-EB074C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified value: 1,03 ± 0,07 mg/kg	Silver (Ag)	[7440-22-4]	0,000103	%		
Certified value: 1,23 ± 0,08 mg/kg	Arsenic (As)	[7440-38-2]	0,000123	%		
Certified value: 0,52 ± 0,06 mg/kg	Gold (Au)	[7440-57-5]	0,000052	%		
Certified value: 0,31 ± 0,06 mg/kg	Beryllium (Be)	[7440-41-7]	0,000031	%		
Certified value: 0,51 ± 0,04 mg/kg	Bismuth (Bi)	[7440-69-9]	0,000051	%		
Certified value: 0,4 ± 0,04 mg/kg	Cadmium (Cd)	[7440-43-9]	0,00004	%		
Certified value: 0,83 ± 0,06 mg/kg	Cobalt (Co)	[7440-48-4]	0,000083	%		
Certified value: 0,37 ± 0,04 mg/kg	Chromium (Cr)	[7440-47-3]	0,000037	%		
Certified value: 5,8 ± 0,8 mg/kg	Iron (Fe)	[7439-89-6]	0,00058	%		
Certified value: 0,49 ± 0,07 mg/kg	Indium (In)	[7440-74-6]	0,000049	%		
Certified value: 2,03 ± 0,27 mg/kg	Magnesium (Mg)	[7439-95-4]	0,000203	%		
Certified value: 0,93 ± 0,07 mg/kg	Manganese (Mn)	[7439-96-5]	0,000093	%		
Certified value: 0,61 ± 0,08 mg/kg	Nickel (Ni)	[7440-02-0]	0,000061	%		
Certified value: 1,53 ± 0,25 mg/kg	Phosphorus (P)	[7723-14-0]	0,000153	%		
Certified value: 2,7 ± 0,4 mg/kg	Lead (Pb)	[7439-92-1]	0,00027	%		
Certified value: 0,57 ±	Antimony (Sb)	[7440-36-0]	0,000057	%		

0,04 mg/kg				
Certified value: 0,55 ±	Selenium (Se)	[7782-49-2]	0,000055	%
0,07 mg/kg				
Certified value: 0,5 ±	Tellurium (Te)	[13494-80-9]	0,00005	%
0,06 mg/kg				
Certified value: 0,97 ±	Titanium (Ti)	[7440-32-6]	0,000097	%
0,18 mg/kg				
Certified value: 2,2 ±	Zinc (Zn)	[7440-66-6]	0,00022	%
0,4 mg/kg				
Indicative value	Mercury (Hg)	[7439-97-6]	<0,00001	%
Indicative value	Sulfur (S)	[7704-34-9]	0,00033	%
Indicative value	Tin (Sn)	[7440-31-5]	0,00015	%
Indicative value	Tungsten (W)	[7440-33-7]	<0,000025	%
Indicative value	Zirconium (Zr)	[7440-67-7]	0,00088	%