

**Single-Element Standard for ICP/MS - Beryllium (Be) 1000 mg/L in 2-5% HNO<sub>3</sub>**

Art. ID LM24-M207.2NP,L1  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	1000	mg/L		Be4O(C2H3O 2)6, 99.99 9%