

**Polyprenol mix C50-C70, Purity: >98%, solution in Hexane**

Art. ID LA-89-5150  
Unit 5 mg  
Deliverydetails De Minimis UN1208 II E2

| Text/Information | Analyte/Parameter | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
|                  | Polyprenol        |              |                     |      |        |        |
|                  | Decaprenol        | [51743-70-5] |                     |      |        |        |
|                  | Undecaprenol      | [50426-50-1] |                     |      |        |        |
|                  | Dodecaprenol      | [10589-57-8] |                     |      |        |        |
|                  | Tridecaprenol     | [52050-19-8] |                     |      |        |        |
|                  | Tetradecaprenol   |              |                     |      |        |        |