

Custom Solution - 32 components in in HNO3 5%, HCl 1% /tr. HF

Art. ID LM24-CUS-70254
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

CRM, ISO 17034 and ISO 17025, Conservation 12 months. (6 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Gold (Au)	[7440-57-5]	10	µg/mL		
	Boron (B)	[7440-42-8]	10	µg/mL		
	Beryllium (Be)	[7440-41-7]	10	µg/mL		
	Cobalt (Co)	[7440-48-4]	10	µg/mL		
	Chromium (Cr)	[7440-47-3]	10	µg/mL		
	Copper (Cu)	[7440-50-8]	10	µg/mL		
	Iron (Fe)	[7439-89-6]	10	µg/mL		
	Germanium (Ge)	[7440-56-4]	10	µg/mL		
	Hafnium (Hf)	[7440-58-6]	10	µg/mL		
	Iridium (Ir)	[7439-88-5]	10	µg/mL		
	Potassium (K)	[7440-09-7]	10	µg/mL		
	Lithium (Li)	[7439-93-2]	10	µg/mL		
	Manganese (Mn)	[7439-96-5]	10	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	10	µg/mL		
	Niobium (Nb)	[7440-03-1]	10	µg/mL		
	Nickel (Ni)	[7440-02-0]	10	µg/mL		
	Osmium (Os)	[7440-04-2]	10	µg/mL		
	Palladium (Pd)	[7440-05-3]	10	µg/mL		
	Platinum (Pt)	[7440-06-4]	10	µg/mL		
	Rhenium (Re)	[7440-15-5]	10	µg/mL		
	Rhodium (Rh)	[7440-16-6]	10	µg/mL		
	Ruthenium (Ru)	[7440-18-8]	10	µg/mL		
	Antimony (Sb)	[7440-36-0]	10	µg/mL		
	Silicon (Si)	[7440-21-3]	10	µg/mL		
	Tin (Sn)	[7440-31-5]	10	µg/mL		
	Tantalum (Ta)	[7440-25-7]	10	µg/mL		
	Tellurium (Te)	[13494-80-9]	10	µg/mL		
	Titanium (Ti)	[7440-32-6]	10	µg/mL		
	Vanadium (V)	[7440-62-2]	10	µg/mL		

Tungsten (W)	[7440-33-7]	10	µg/mL
Zinc (Zn)	[7440-66-6]	10	µg/mL
Zirconium (Zr)	[7440-67-7]	10	µg/mL