

ARCONIC - Aluminium base, RM, Alloy: 2199, cast, Batch A, Disc 64 mm Ø x 25 mm

Art. ID ALCOA-Li-2199
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 0,0001 | % | | |
| | Arsenic (As) | [7440-38-2] | ~0,0001 | % | | |
| | Boron (B) | [7440-42-8] | 0,0001 | % | | |
| | Calcium (Ca) | [7440-70-2] | 0,0003 | % | | |
| | Chromium (Cr) | [7440-47-3] | 0,0006 | % | | |
| | Copper (Cu) | [7440-50-8] | 2,58 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,0352 | % | | |
| | Gallium (Ga) | [7440-55-3] | 0,005 | % | | |
| | Lithium (Li) | [7439-93-2] | 1,63 | % | | |
| | Magnesium (Mg) | [7439-95-4] | 0,201 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,301 | % | | |
| | Sodium (Na) | [7440-23-5] | 0,0002 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,0004 | % | | |
| | Phosphorus (P) | [7723-14-0] | ~0,0001 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,0005 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,023 | % | | |
| | Strontium (Sr) | [7440-24-6] | 0,0001 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,048 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,0064 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,55 | % | | |
| | Zirconium (Zr) | [7440-67-7] | 0,0837 | % | | |