

**ARCONIC - Aluminium base, RM, Alloy: 2199, cast, Batch A, Disc 64 mm Ø x 25 mm**

Art. ID ALCOA-Li-2199  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0001	%		
	Arsenic (As)	[7440-38-2]	~0,0001	%		
	Boron (B)	[7440-42-8]	0,0001	%		
	Calcium (Ca)	[7440-70-2]	0,0003	%		
	Chromium (Cr)	[7440-47-3]	0,0006	%		
	Copper (Cu)	[7440-50-8]	2,58	%		
	Iron (Fe)	[7439-89-6]	0,0352	%		
	Gallium (Ga)	[7440-55-3]	0,005	%		
	Lithium (Li)	[7439-93-2]	1,63	%		
	Magnesium (Mg)	[7439-95-4]	0,201	%		
	Manganese (Mn)	[7439-96-5]	0,301	%		
	Sodium (Na)	[7440-23-5]	0,0002	%		
	Nickel (Ni)	[7440-02-0]	0,0004	%		
	Phosphorus (P)	[7723-14-0]	~0,0001	%		
	Lead (Pb)	[7439-92-1]	0,0005	%		
	Silicon (Si)	[7440-21-3]	0,023	%		
	Strontium (Sr)	[7440-24-6]	0,0001	%		
	Titanium (Ti)	[7440-32-6]	0,048	%		
	Vanadium (V)	[7440-62-2]	0,0064	%		
	Zinc (Zn)	[7440-66-6]	0,55	%		
	Zirconium (Zr)	[7440-67-7]	0,0837	%		