

ARCONIC - Aluminium base, RM, Alloy: 7085, cast, typical analysis, ppm values not shown, Disc 64 mm Ø x 25 mm

Art. ID ALCOA-SS-7085
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chromium (Cr)	[7440-47-3]	0,02	%		
	Copper (Cu)	[7440-50-8]	1,7	%		
	Iron (Fe)	[7439-89-6]	0,04	%		
	Magnesium (Mg)	[7439-95-4]	1,5	%		
	Manganese (Mn)	[7439-96-5]	0,02	%		
	Silicon (Si)	[7440-21-3]	0,03	%		
	Titanium (Ti)	[7440-32-6]	0,03	%		
	Zinc (Zn)	[7440-66-6]	7,5	%		
	Zirconium (Zr)	[7440-67-7]	0,1	%		