

ARCONIC - Aluminium base, RM, Alloy: 4001, cast, Batch A, Disc 64 mm Ø x 25 mm

Art. ID ALCOA-SS-SI4001
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	~0,0005	%		
	Boron (B)	[7440-42-8]	~0,002	%		
	Bismuth (Bi)	[7440-69-9]	~0,0003	%		
	Calcium (Ca)	[7440-70-2]	~0,02	%		
	Cadmium (Cd)	[7440-43-9]	~0,0001	%		
	Cobalt (Co)	[7440-48-4]	~0,0008	%		
	Chromium (Cr)	[7440-47-3]	0,0046	%		
	Copper (Cu)	[7440-50-8]	0,012	%		
	Iron (Fe)	[7439-89-6]	0,33	%		
	Gallium (Ga)	[7440-55-3]	0,008	%		
	Lithium (Li)	[7439-93-2]	~0,0003	%		
	Magnesium (Mg)	[7439-95-4]	0,054	%		
	Manganese (Mn)	[7439-96-5]	0,014	%		
	Sodium (Na)	[7440-23-5]	~0,0001	%		
	Nickel (Ni)	[7440-02-0]	0,0085	%		
	Phosphorus (P)	[7723-14-0]	~0,012	%		
	Lead (Pb)	[7439-92-1]	~0,001	%		
	Scandium (Sc)	[7440-20-2]	~0,006	%		
	Silicon (Si)	[7440-21-3]	38,1	%		
	Tin (Sn)	[7440-31-5]	~0,0002	%		
	Strontium (Sr)	[7440-24-6]	~0,007	%		
	Titanium (Ti)	[7440-32-6]	0,02	%		
	Vanadium (V)	[7440-62-2]	0,01	%		
	Zinc (Zn)	[7440-66-6]	0,012	%		
	Zirconium (Zr)	[7440-67-7]	~0,005	%		