

## **Monomethylarsonic Acid Standard Solution**

Art. ID	NIST-3030
Unit	2 x 5 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM) is intended for use as a primary calibration standard for the quantitative determination of monomethylarsonic acid (MMA). A unit of NIST-3030 consists of two 5 mL sealed borosilicate glass amber ampoules of water solution prepared gravimetrically to contain a known mass fraction of MMA. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified Mass Fraction	Monomethylarsonic acid	[124-58-3]	17,64 ± 0,15	mg/kg		
Value of MMA (as As)	(MMA)					