

Dimethylarsinic Acid Standard Solution

Art. ID	NIST-3031
Unit	2 x 5 mL
Deliverydetails	No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM) is intended for use as a primary calibration standard for the quantitative determination of dimethylarsinic acid (DMA). A unit of NIST-3031 consists of two 5 mL sealed borosilicate glass amber ampoules of water solution prepared gravimetrically to contain a known mass fraction of DMA. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified Mass Fraction	Dimethylarsinic acid (D		20,47 ± 0,18	mg/kg		
Value of DMA (as As)	MA)					