

Custom Solution - 9 components in HNO₃ (2%) (CRM, ISO 17034 and ISO 17025)

Art. ID LM24-CUS-70259
Unit 200 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	10	mg/L		
	Arsenic (As)	[7440-38-2]	1	mg/L		
	Iron (Fe)	[7439-89-6]	20	mg/L		
	Lead (Pb)	[7439-92-1]	10	mg/L		
	Copper (Cu)	[7440-50-8]	10	mg/L		
	Zinc (Zn)	[7440-66-6]	10	mg/L		
	Manganese (Mn)	[7439-96-5]	10	mg/L		
	Cobalt (Co)	[7440-48-4]	10	mg/L		
	Nickel (Ni)	[7440-02-0]	10	mg/L		