

## **Arsenobetaine Standard Solution**

Art. ID	NIST-3033
Unit	2 x 5 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM) is intended for use as a primary calibration standard for the quantitative determination of arsenobetaine (AB). A unit of NIST-3033 consists of two 5 mL sealed borosilicate glass amber ampoules of water solution prepared gravimetrically to contain a known mass fraction of AB. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified Mass Fraction	Arsenobetaine (AB)	[64436-13-1]	19,06 ± 0,27	mg/kg		
Value of AB (as As)						