

Arsenocholine Standard Solution

Art. ID	NIST-3034
Unit	2 x 5 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Description

This Standard Reference Material (SRM) is intended for use as a primary calibration standard for the quantitative determination of arsenocholine (AC). A unit of NIST-3034 consists of two 5 mL sealed borosilicate glass amber ampoules of water solution prepared gravimetrically to contain a known mass fraction of AC. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified Mass Fraction	Arsenocholine (AC)		19,77 ± 0,27	mg/kg		
Value of AC (as As)						