

lipase Lipex 100L (Lipex 100L1), covalent on IB-150A. Activity: 3000 TBU / g, 9000 PLU / g

Art. ID	CVI-IMMLIPX-T2-150-100G
Unit	100 g
Deliverydetails	No Dangerous Good

Description

Unit definition: 1 TBU unit = 1 μ mol butyric acid released per minute / g immobilized enzyme at 40°C and pH 7.5 ///

Unit definition: 1 PLU unit = 1 μ mol propyl laurate formed per minute / g immobilized enzyme at 60°C and pH 7.5.