

lipase from Aspergillus niger (Lipase A4), covalent on IB-150A. Activity: 40 TBU / g

Art. ID	CVI-IMMANA-T2-150-1KG
Unit	1 kg
Deliverydetails	No Dangerous Good

Description

Unit definition: 1 TBU unit = 1 μ mol butyric acid released per minute / g immobilized enzyme at 40°C and pH 7.5 ///

Unit definition: 1 PLU unit = 1 μ mol propyl laurate formed per minute / g immobilized enzyme at 60°C and pH 7.5.