

**lipase Resinase HT (Resinase HT1), covalent on IB-150A. Activity: 3000 TBU / g, 6000 PLU / g**

Art. ID	CVI-IMMRES-T2-150-1KG
Unit	1 kg
Deliverydetails	No Dangerous Good

Description

Unit definition: 1 TBU unit = 1  $\mu$ mol butyric acid released per minute / g immobilized enzyme at 40°C and pH 7.5 ///

Unit definition: 1 PLU unit = 1  $\mu$ mol propyl laurate formed per minute / g immobilized enzyme at 60°C and pH 7.5.