

**lipase Lipex 100L (Lipex 100L1), covalent on IB-150A. Activity: 3000 TBU / g, 9000 PLU / g**

|                 |                         |
|-----------------|-------------------------|
| Art. ID         | CVI-IMMLIPX-T2-150-10KG |
| Unit            | 10 kg                   |
| Deliverydetails | No Dangerous Good       |

**Description**

Unit definition: 1 TBU unit = 1  $\mu$ mol butyric acid released per minute / g immobilized enzyme at 40°C and pH 7.5 ///  
Unit definition: 1 PLU unit = 1  $\mu$ mol propyl laurate formed per minute / g immobilized enzyme at 60°C and pH 7.5.