

lipase from Bacillus subtilis (lipase7) adsorbed on dry acrylic beads Activity: 2000 TBU / g

Art. ID	CVI-IMMAULI-T1-350-10KG
Unit	10 kg
Deliverydetails	No Dangerous Good

Description

Unit definition: 1 TBU unit = 1 μ mol butyric acid released per minute / g immobilized enzyme at 40°C and pH 7.5 ///
Unit definition: 1 PLU unit = 1 μ mol propyl laurate formed per minute / g immobilized enzyme at 60°C and pH 7.5.