

## **Sulfur and Btu In Paraffin Oil Reference Material**

Art. ID	AR-000
Unit	0,1 L
Deliverydetails	No Dangerous Good

### Description

AR-000 is a paraffin oil certified reference material (CRM). The intended use of this CRM is for the verification and calibration of instrumental methods for the determination of sulfur and heat of combustion (see references below). This CRM can also be used to validate value assignment of in-house reference materials. A unit consists of one bottle containing 100 mL of reference material as a liquid. All reference materials should be verified as fit for purpose prior to use.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]	0,0005 ± 0,0002	%		
	Calorific Value		19887	BTU/lb		