

Custom Solution - 3 components in 5% HNO₃ (ISO 17034 / ISO 17025)

Art. ID LM24-CUS-70269
Unit 50 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Stored under normal laboratory conditions. Conservation 24 months.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	750	mg/L		
	Thallium (Tl)	[7440-28-0]	400	mg/L		
	Lead (Pb)	[7439-92-1]	250	mg/L		