

**Tiegel (Crucible, with pins at reinforced ring, TG2, Pt/Au 95/5)**

Art. ID XRF-SC025A044RF  
Unit approx. 45 g  
Deliverydetails No Dangerous Good

**Description**

\*\*Used on the original Schoeps fusion machines as well as on machines of type HD-VAA, HD-Vulcan and Autofluxer\*\*

FAB - Crucible Pt95/Au5 25ml SC6 3 Pin with ring, reinforced with spout

RFA Tiegel aus PtAu5

passend für Autofluxer, Schoeps, Vulcan

mit verstärktem Rand und Halterung

Innendurchmesser oben 34,00 mm

Innendurchmesser unten 20,00 mm

Höhe 38,00 mm

Bodendicke 0,50 mm

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Platinum (Pt)	[7440-06-4]	95	%		
	Gold (Au)	[7440-57-5]	5	%		