

FAB - Casting Dish XRS Pt95/Au5 39/41x3mm SM5 / RFA Abgiessschale aus Platin/ Gold 95/5 - XRS

Art. ID XRF-SM3941AX045DHF
Unit approx. 45 g
Deliverydetails No Dangerous Good

Description

Used on the original Schoeps fusion machines as well as on machines of type HD-VAA, HD-Vulcan and Autofluxer

FAB - Casting Dish XRS Pt95/Au5 39/41x3mm SM5

RFA Abgiessschale aus Platin/ Gold 95/5 - XRS

passend für Autofluxer, Schoeps, Vulkan

Durchmesser 39/ 41 mm

Kantenlänge 51,00 mm

Innenhöhe 3,00 mm

Bodendicke 0,80 mm

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Platinum (Pt)	[7440-06-4]	95	%		
	Gold (Au)	[7440-57-5]	5	%		