

FAB - Lid Pt95/Au5 to Suit 20ml Std Crucible GCS20L Deckel aus Platin/ Gold 95/5

Art. ID	XRF-GLS020A004LF
Unit	approx. 4,7 g
Deliverydetails	No Dangerous Good

Description

Used on the original Schoeps fusion machines as well as on machines of type HD-VAA, HD-Vulcan and Autofluxer

FAB - Lid Pt95/Au5 to Suit 20ml Std Crucible GCS20L

Deckel aus Platin/ Gold 95/5

passend für Tiegel SC6

Durchmesser 33,00 mm

Blechdicke 0,15 mm

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Platinum (Pt)	[7440-06-4]	95	%		
	Gold (Au)	[7440-57-5]	5	%		