

## **Silver nitrate - Assay/isotopic composition**

Art. ID	NIST-978a
Unit	0,25 g
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1479 II E2

### Description

Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Purity	Silver nitrate	[7761-88-8]	99,99	%		
Abundance ratio	107Ag/109Ag		1,07638			
Isotopic composition	107Ag		51,8392	%		
Isotopic composition	109Ag		48,1608	%		