

N-Acetyl-D-alanyl-D-alanine

Art. ID TRC-A168068-250MG
Unit 250 mg
Deliverydetails De Minimis

Description

Application Note: N- γ -Acetyl- γ -D- γ -alanyl-D-alanine was used as the target in the target-accelerated combinatorial synthesis of vancomycin dimers. D-Ala-D-Ala possess the greatest affinity for vancomycin. Vancomycin and N- γ -Acetyl- γ -D- γ -alanyl-D-alanine are used as a model system to demonstrate the use of affinity capillary eletrophoresis.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N- γ -Acetyl- γ -D- γ -alanyl-D-alanine	[19993-26-1]				