

## **2-Allyl-2-phenyl-4-pentenitrile**

Art. ID	TRC-A549215-2.5G
Unit	2,5 g
Deliverydetails	De Minimis

### Description

Application Note: 2-Allyl-2-phenyl-4-pentenitrile is an intermediate in synthesizing alpha-Phenyloxiranepropanenitrile (P335835), used in the preparation of CC-1065 and its analogs.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Allyl-2-phenyl-4-pent enenitrile	[28049-67-4]				