

Angiotensin I (human)

Art. ID NIST-998
Unit 0,5 mg
Deliverydetails No Dangerous Good (Dry ice shipment) /not restricted

Description

This Standard Reference Material (SRM®) is certified as a chemical of known purity. It is intended primarily for use in the calibration and standardization of the renin assay and as a reference peptide for amino acid analysis and high performance liquid chromatography (HPLC). Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Angiotensin I		94,1 ± 0,9	%		