

(2S)-2-Amino-4-azidobutanoic Acid Hydrochloride

Art. ID TRC-A584003-50MG
Unit 50 mg
Deliverydetails De Minimis

Description

Application Note: (2S)-2-Amino-4-azidobutanoic Acid Hydrochloride is a modified amino acid used for the introduction of bioorthogonal handles into proteins. It is used in the synthetic preparation of glucosyl phosphine derivatives for site-selective traceless Staudinger ligation with azide containing proteins.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(2S)-2-Amino-4-azidobutanoic Acid Hydrochloride	[942518-29-8]				