

## **Cigarette ignition strength standard**

Art. ID	NIST-1082
Unit	1 carton (200 cigarettes)
Deliverydetails	Dry ice shipment /not restricted

### **Description**

This Standard Reference Material (SRM®) is intended for use by test laboratories to assess and control their testing of cigarette ignition strength in accordance with ASTM Standard Methods E 2187-04 (or ASTM E2187-02b). The SRM® unit consists of one carton of cigarettes containing 10 packs of 20 cigarettes each. Certified value /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ignition Strength		12.6 % ± 3.3 %		ASTM E 218 7-04	