

Amoxicillin Dimer tri-Sodium Salt (penicilloic acid form) (Mixture of Diastereomers)

Art. ID	TRC-A634520-0.5MG
Unit	0,5 mg
Deliverydetails	De Minimis

Description

Application Note: Amoxicillin dimer tri-Sodium Salt (penicilloic acid form) is an impurity of Amoxicillin (A634235), a beta-lactam family antibiotic that is commonly used in conjunction with clavulanic acid (C56370) to the treat bacterial infections.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Amoxicillin Dimer tri-Sodium Salt (penicilloic acid form) (Mixture of Diastereomers)					