

(D-Asn)-Desmopressin

Art. ID TRC-A788858-5MG
Unit 5 mg
Deliverydetails De Minimis

Description

Application Note: (D-Asn)-Desmopressin is an impurity of Desmopressin (D296790), a synthetic analog of Vasopressin with antidiuretic activity and decreased pressor effects. Improves hemostasis by increasing plasma levels of von Willebrand factor, factor VIII and tissue plasminogen activator. Antidiuretic, hemostatic.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(D-Asn)-Desmopressin					