

3-Azido-7-hydroxycoumarin

Art. ID TRC-A850625-5MG
Unit 5 mg
Deliverydetails De Minimis

Description

Application Note: 3-Azido-7-hydroxycoumarin, is a Fluorescent dye that has shown to quickly and efficiently label proteins in bicarbonate buffered solns. (Abs/Em = 404/477 nm)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-Azido-7-hydroxycoumarin	[817638-68-9]				
	in					