

## **Aluminum Alloy 6011**

Art. ID	NIST-1258-I
Unit	disc
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM®) is intended primarily for use in optical emission and X-ray spectrometric methods of analysis. A unit of NIST-1258-I consists of wrought and annealed disks, 35 mm (1 3/8 in) diameter and 19 mm (3/4 in) thick. Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	0,8	%		
	Iron (Fe)	[7439-89-6]	0,08	%		
	Copper (Cu)	[7440-50-8]	0,848	%		
	Manganese (Mn)	[7439-96-5]	0,481	%		
	Magnesium (Mg)	[7439-95-4]	1	%		
	Zinc (Zn)	[7440-66-6]	1,03	%		