

## **Fibrous glass blanket for high precision measurements - Thermal resistance**

Art. ID	NIST-1452
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM®) is intended for use in evaluation of a guarded hot plate (GHP) or the calibration of a heat flow meter (HFM). It is supplied as a fibrous glass batt of nominal dimensions 60 x 60 x 2.54 cm. Each unit of NIST-1452 is a individually characterised specimen.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Please aks for the current certificate						