

Out of Stock - Item is not available at this time - Kerosene - Sulfur

Art. ID	NIST-1616b
Unit	100 mL
Deliverydetails	No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM®) is intended for use in the evaluation of methods and the calibration of instruments used in the determination of total sulfur in fuel oils or materials of a similar matrix. NIST-1616b consists of a special low sulfur kerosene (No. 1-K) suitable for use in non flue-connected burner appliances and for use in wick-fed illuminating lamps, as described in ASTM D 3699-92 Specification for Kerosene. A unit of NIST-1616b consists of 100 mL of the low sulfur kerosene in an amber glass bottle. Certified value /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]	8,41	mg/kg		