

Anode tin - Lead, Block

Art. ID	NIST-1727
Unit	disc
Deliverydetails	No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM®) is intended primarily for use in evaluating chemical methods and instrumental methods of analysis of refined pig tin in anode form as specified by ASTM International B 339-00 Standard Specification for Pig Tin. A unit of NIST-1727 consists of a block of Grade ??A?? tin for the Manufacture of Tinplate with dimensions approximately 30 mm x 30 mm x 30 mm. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lead (Pb)	[7439-92-1]	0,003326	%		