

## **Zinc-Aluminum Alloy**

Art. ID	NIST-1741
Unit	disc
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM®) is in the form of a disk, approximately 50.8 mm (2 in) in diameter and 12.7 mm (0.5 in) thick intended for use with optical emission and X-ray fluorescence spectrometric methods of analysis. This material is one in a series of seven zinc base alloys prepared to cover a range of aluminum and lead compositions of interest to the zinc and galvanizing industries. The other zinc alloys in this series are SRMs 1736, 1737, 1738, 1739, 1740, and 1742, all in disk form, SRM 2139 is the same material as SRM 1739 but supplied in chip form. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,5242	%		
	Lead (Pb)	[7439-92-1]	0,1571	%		