

## **Linerboard - Tape adhesion testing (sheet form)**

Art. ID	NIST-1810a
Unit	50 each
Deliverydetails	No Dangerous Good /not restricted

### **Description**

This Standard Reference Material (SRM®) is intended to provide a uniform source of linerboard for use with ASTM D 2860 Standard Test Method for Adhesion of pressure-sensitive tape to fiberboard at 90 degree angle and constant stress. Each unit consists of fifty 21.6 cm x 28 cm (8.5 in x 11 in) linerboard sheets of Mosinee 696-C paper, sealed in a moisture resistant foil-lined package.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Please ask for current certificate						