

Lubricating base oil - Chlorine

Art. ID	NIST-1818a
Unit	set (5)
Deliverydetails	No Dangerous Good /not restricted

Description

Set of 5 x 20 mL ampoules This Standard Reference Material (SRM) is intended for use as an analytical standard for the determination of total chlorine in lubricating base oils or materials of similar composition. NIST-1818a consists of a series of five 20 mL ampoules of lubricating base oils at different chlorine concentrations. Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Level 1818a- I	Chlorine (Cl)	[7782-50-5]	31,6	mg/kg		
Level 1818a-II	Chlorine (Cl)	[7782-50-5]	60	mg/kg		
Level1818a-III	Chlorine (Cl)	[7782-50-5]	78,2	mg/kg		
Level1818a-IV	Chlorine (Cl)	[7782-50-5]	154,4	mg/kg		
Level1818a-V	Chlorine (Cl)	[7782-50-5]	234	mg/kg		