

Silicon nitride - Specific surface area

Art. ID	NIST-1900
Unit	4 g
Deliverydetails	No Dangerous Good /not restricted

Description

Intended for use in the calibration of BET instruments used to measure specific surface area (SSA) in the range 0.1 m²/g to 1000 m²/g. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Specific surface area (BET) multipoint		2,85	m ² /g		
	Specific surface area (BET) single point		2,79	m ² /g		