

## **Polystyrene (0.5 wt. % in water) - Particle size**

Art. ID	NIST-1961
Unit	5 mL
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM®) is intended for use as a primary particle size reference standard for the calibration of particle measuring instruments including flow-through counters and optical and electron microscopes. NIST-1961 is a suspension of polystyrene spheres in water at a weight concentration of about 0,5%. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Number average diameter		29,64	µm		