

Out of Stock - Product is not available at this time - Low-Energy Charpy V-Notch Specimens (Self-Verification, 8-mm Striker) (ASTM E23)

Art. ID	NIST-2093
Unit	set (5)
Deliverydetails	No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM) is intended for the verification of Charpy machines equipped with an 8-mm striker, in accordance with the current ASTM Standard E23 and the current ISO Standard 148-2. Each SRM consists of a set of five specimens needed to perform one verification. NIST-2093 is used for in-house verification (self-verification).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Please ask for current certificate						