

## **LA Steel, Carbon & Sulfur only**

Art. ID	NIST-2160
Unit	200 g
Deliverydetails	No Dangerous Good /not restricted

### Description

This Standard Reference Material (SRM) is in the form of pins that are 4 mm in diameter, 12 mm long and weigh approximately 1g each. SRM 2160 is primarily intended for use in calibrating instruments used in the determinations of carbon und sulfur in steel. Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,584	%		
	Sulfur (S)	[7704-34-9]	0,012	%		