

## **High purity Iron (chips)**

Art. ID                      NIST-2168  
Unit                        150 g  
Deliverydetails            No Dangerous Good /not restricted

### **Description**

This Standard Reference Material (SRM®) is high-purity iron in the form of chips and is intended primarily for evaluation of methods for analysis of trace elements in iron and materials of similar matrix. It should not be used for calibration. A unit of NIST-2168 consists of one bottle containing approximately 150 g of chips. Certified values /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]	10,7 ± 2,8	mg/kg		